### **Attachment**

## YC28 VENTILATORS



- Built in oxygen generator, high flow output, and no need for external oxygen supply
- Can effectively work on acute respiratory failure symptoms (ARDS) and hypoxic symptoms.
- Humidification effect can reach h80%.
- One stop structure, easy to move and operate.

#### Model No. YC28 Key Specifiation:

- 1.Pressure rangeL 4-30cm H20
- 2.Treatment modes: CPAPS.S-Auto,TS/T
- 3.Flow range: 20-80L/min
- 4.Temperature range: 30C-37C
- 5.Oxygen flow: 21-66%6.Humidification: 50-80%
- 7. Parameter display: Oxygen density,

oxygen flow, output flow

#### Application:

Pneumology, Emergency ICU, Neurosurgery, Neurology, Sleep clinics, Chest Pain Center etc.

THE PROPERTY AND ADDRESS OF THE PROPERTY AND IN COLUMN THE PROPERTY AND ADDRESS OF THE PROPERTY ADDRESS OF THE









此报价有效期至2020-3-29 Special Quotation from Yamind,valid to 2020-3-29.										
型号 ITEM	产品图片 Product picture	配件: Accessories:	产品配置 Configuration	单价(未税) Price 500-1000PCS						
YC-28 正圧通气呼吸 制氧一体系统		配件: Accessories: 認化器: 管路、电源、全验面罩、适配器、4Gb SD卡、说明书、电源线、操作指弯、装箱单、便预包、过途根、导管等。实化接头、Humddiffe.Af Tobing Air Filler, Power Supplymask M,Adaptor,4Gb SD Card, User Manual,Carrying bag Power Cable,Quick Start Guide,Packing List,Nasal oxygen tube.Atomizing joint.	DM28-30ST-B机器一会.10L制氧机一合。治疗模式CPAP, APAP, S, S-Auto, T, S/T.AVAPS模式: 压力范围-30 CmH20/劳女计数模算/发动设计疗曲线量元/0.2cmH20压力精度/投管功能/WIFI诗数据连接/WIFI46B SD数据存储/标准包装配件* DM28-30ST-B*1ps:70LOxygen concentrator*1pcs.Treatment modes: CPAP, APAP, S, S-Auto, T, S/T, AVAPS control/Pressure range:4-30 cmH2O/5* LCD touch screen/Real time curve/0.2cmH2O pressure increment/alarm functions/WIFI/4GB SD card/standard package accessories  DM28T: *Size: 380*170*335mm *Weight: 3.9KG CPS01: *Size: 475*375*615mm *Weight: 16.5KG							

保修期:主机保修期2年,加湿器保修期1年,加湿器保修期3个月用于管道、面罩、水箱等附件。

Warranty Period: 2-year warranty for the host, 1-year warranty for the humidifier and 3-month warranty

for accessories such as the tubing, mask, and water tank. 备注:1.图片仅供参考,

付款方式: 100%电汇预付。

Payment Terms: 100% T/T advanced.

交货时间: 收到全部贷款后,100台以内,15天内交付,大于100台,需排期15-45天内(疫情期间视情况而定)。装运条件EXW,可通过快递、空运或海运到岸价指定。提醒:产品图片仅供参考。 Lead Time: Within 100 sets, delivered within 7 days, more than 100 sets Scheduled period 15-45 days after receipt of full purchase payment(Lead time need to confirm during the epidemic situation

period).
Shipment Terms: EXW, and could be appointed on CIF basis by courier, air or sea. Remind: Product pictures are for reference only. 深圳市大雅医疗技术有限公司

© 2019 Ymind All rights reserved.ShenZhen Yamind Medical Technology Co., Ltd.




## Shenzhen Yamind medical


## YAMIND BiPAP & Oxygen Concentrator Combo Packlist 装箱单/重量单


Sold to	Messrs:	No.:								
		Invoice No.:								
				Date:						
		L/C No.:								
				S/C No.:						
From		to			by					
Part- list	Part-Name	Quantity		Note 备注	Packaging					
1	部件名称 DM28-30ST-B BiPAP Host	<u>数量</u> 1	Vo	ntilator Host						
2	Air Tubing	1	•	initiator riost						
3		2	<del>                                     </del>		Packaging-1					
_	Air Filter				机器净重(KG)Net weight:					
4	Power Cable	1			*Size: 271.5*143*139.5mm *Weight: 1.8KG					
5	Adaptor	1			包装毛重The packing gross weight:					
6	mask M	1	1		*Size: 380*170*335mm					
7	4Gb SD Card	1			*Weight: 3.9KG 五合一尺寸(mm)Five-in-one size: *Size: 930*405*365mm					
8	User Manual	1								
9	Carrying bag	1			五合一毛重Five-in-one gross weight: *Weight: 21KG					
10	Quick Start Guide	1			Weight. 21KG					
11	Packing List	1								
12	Oxygen Concentrator host	1	1 Oxygen		Packaging-1 机器净重(KG)Net weight:					
40	0.57.50	1	Con	centrator hos	*Size: 400*300*540					
	<ul> <li>Oxygen Concentrator Air Filter</li> <li>Nasal oxygen tube</li> </ul>				*Weight: 14.5 包装毛重The packing gross weight:					
14										
15	Atomizing joint	1			*Size: 475*375*615mm					






# CPAP & Oxygen Concentrator Combo





1. Plug the CPAP power cable and PS/2 data cable 2. Connect T- connector with the CPAP outlet. into the power interface and PS/2 port.



3. Set the catheter over the oxygen outlet tightly, then 4. Connect the air tube to the T- connector set the other end of the catheter to the CPAP outlet.